

How to insert a battery

Step 1

Open battery door by pushing lip down.

Step 2

When you open the battery door fully it will look like this.

Step 3

Place the battery into the slot with the orange tab facing up.



Step 4

Peel the orange tab off the battery.

Step 5

To switch the hearing aid on, push the battery door all the way in so it clicks shut.



Leaflet No: 01-14-101 12/2014 Review date: 12/2016



Yeovil District Hospital MHS

NHS Foundation Trust





How to insert a battery

Step 1

Open battery door by pushing lip down.

Step 2

When you open the battery door fully it will look like this.

Step 3

Place the battery into the slot with the orange tab facing up.

Step 4

Peel the orange tab off the battery.

Step 5

To switch the hearing aid on, push the battery door all the way in so it clicks shut.

Leaflet No: 01-14-101 12/2014 Review date: 12/2016













